

Microkleen™

Aldehyde-Free, Quat-Based Rapid Disinfectant

Composition :

Active ingredients :

Propanol-1	30 g
Propanol-2	45 g
Didecyl dimethyl ammonium chloride	0.5 g
N, alkyl dimethyl benzyl ammonium chloride	0.5 g
Poly-hexamethylene biguanide hydro-chloride	0.5 g
Corrosion inhibitors, Chelating agents	

About Microkleen

Microkleen is a ready-to-use, disinfectant cleaner and deodorant designed for rapid disinfecting and controlling mold and mildew on hard, nonporous surfaces. Use this non-abrasive product to quickly remove dirt, grime, mold, mildew and other common soils found in hospitals, nursing homes, clinics, medical offices, laboratories, veterinary clinics, food processing plants and pharmaceutical industries.

Quat-based, hospital-grade disinfectant cleaner that cleans, disinfects and deodorizes in one step.

Broad Spectrum Activity

- ▶ Highly effective against a wide variety of pathogenic microorganisms including bacteria, antibiotic resistant bacteria, viruses, fungi, mold and mildew.
- ▶ Bactericidal, Fungicidal, Mildewstatic & Virucidal

Effective

- ▶ Ready-to-use liquid disinfects in 30 seconds to five minutes contacts.
- ▶ Cleans and disinfects surfaces
- ▶ Quickly removes dirt, grime, mold, mildew and other common soils from hard, nonporous surfaces

Versatile

- ▶ Use on hard, nonporous surfaces including vinyl, painted surfaces, plastic, glazed ceramic, glazed porcelain, chrome, laminated surfaces and baked enamel surfaces
- ▶ Effective on stains in health care, educational, lodging and food service environments



Micr**gen**

Redefining Hygiene...



Directions for use :

- ▶ Wipe the surfaces to be disinfected, with a sufficient amount of ready-to-use solution, ensuring complete coverage. The product has got prolonged residual effect. While disinfecting instruments ensure that instruments is wiped free of **Microkleen** with sterile cloth. Wipe the surface to be used with sterile cloth after applying **Microkleen**.

Product specification

Description

Colour Form	Colourless
Fragrance	Odourless
Shelf life	36 months
Solubility	Soluble
Pack Size	100 ml. 500 ml

A product of :

Microgen Hygiene Pvt. Ltd.

313, Midas, 3rd Floor, Sahar Plaza, Andheri-Kurla Road, Andheri (East). Mumbai - 400 059.

Tel : +91 22 2825 1870 **Fax :** +91 22 2832 0335

Website : www.microgenindia.co

WIDE USAGE

For usage in NICU

- ▶ Syringe pump
- ▶ Infusion pump
- ▶ ECG monitor
- ▶ Pulse oxymeter
- ▶ Photo therapy unit
- ▶ Ventilator
- ▶ Accessories – cables, cords, adaptors, connectors, wires
- ▶ Incubators & warmers

For usage in Dental homes

- ▶ Dental chair, equipment & accessories, waste receivers
- ▶ Hand pieces in kidney plates
- ▶ Spitting areas of dental chairs

Usage in operation theatres

- ▶ Operation tables
- ▶ Overhead lights
- ▶ Operation theatres between surgeries
- ▶ Monitors
- ▶ Fibrillator
- ▶ Heart-lung machine

Usage in AKD

- ▶ Haemodialysis system
- ▶ Blood pump
- ▶ Respirator Monitor

Usage in Ophthalmology

- ▶ Tornometers, microscopes, scan probes etc.

Others

- ▶ Imaging machines
- ▶ C.T. scanners
- ▶ HIV endangered areas
- ▶ Pharmaceutical sterile sections