

# D-125<sup>®</sup> Wet Wipes

## Presoaked Surface Disinfection Wipes

### Revolutionary Formula :

- ✓ Aldehyde Free
- ✓ Alcohol Free
- ✓ Bleach Free

### 3rd Generation Quaternary Ammonium Compound

- Disinfectant • Bactericidal
- Virucidal • Fungicidal
- Mildewcidal • Deodorizer

### About :

**D-125 Wet Wipes** are non-woven disposable wipes, presoaked with D-125 solution. D-125 is an US-EPA registered hospital-grade disinfectant. **D-125 Wet Wipes** are intended for use in health care settings such as Hospitals, Speciality health care facilities, Nursing home, Ophthalmic clinics, Dental clinics, Laboratories, Veterinary facilities & many more...


- Universal • Comfortable • Time saving

### The Benefits of Ready to Use Wipes :

The traditional towel and bucket method of environmental cleaning is being replaced in many hospitals by ready to use disinfectant wipes. The pre-soaked disposable wipes are increasingly used for disinfection of near-patient surfaces and prevent the spread of microorganisms and the emergence of nosocomial infections. The report explores the benefits of ready to use wipes as well as reviews the key factor that impact wipes efficacy.

*Reference : Infection Control Today - Special Report August 2015*

**Microgen<sup>®</sup>**  
Redefining Hygiene...



### Features and Benefits :

- Alcohol free, Aldehyde free & Bleach free formula
- Presoaked with D-125 solution, at dilution 1:32
- Cleans and disinfects – one step
- Contact time is 10 minutes or less\*\*
- Shelf life - 90 days (after opening)
- US - EPA registered; meets CDC & OSHA guidelines – Maintain compliance to protocols
- **Diamond-embossed material** – is thick and strong to help clean soiled surfaces

# D-125<sup>®</sup> Wet Wipes

Alcohol, Aldehyde and Bleach Free Formula :

- Bactericidal, fungicidal, and virucidal.
- Cleans and disinfects in one step.
- Tested effective against over 148 types of microorganisms.
- Contact time is 10 minutes or less.
- Use a wipe to remove heavy soil. Unfold a clean wipe & thoroughly wet surface. Treated surface must remain visibly wet for a full 10 minutes. Use additional wipes if needed to assure continuous 10 minutes wet contact time. Let air dry.

## Active Ingredients :

Alkyl Dimethyl Benzyl Ammonium Chloride  
Alkyl Dimethyl Ethylbenzyl Ammonium Chloride  
Inert Ingredients

At dilution 1:32

## Microbiological Efficacy :

**D-125 Wet Wipes Active against Bacteria, Viruses, and Fungi, including :**

- Escherichia coli strain O157:H7 • Pseudomonas aeruginosa • Klebsiella pneumoniae
- Acinetobacter baumannii • Vaccinia Virus
- HIV-1 • Herpes simplex type 2 • Hepatitis B
- Influenza A2-Asian Virus • Candida albicans
- Adenovirus Type 2 • Trichophyton mentagrophytes • Salmonella choleraesuis
- Staphylococcus aureus • Methicillin Resistant Staphylococcus aureus (MRSA) • Vancomycin Resistant Enterococcus (VRE) etc.



**US EPA - Registered Hospital Disinfectant**

(Reg 61178-1) approved as a Hospital Disinfectant with Broad Spectrum Antimicrobial Activity in 10 min.



**WHO Certified**

World Health Organization Approved D-125 for H1N1, H5N1, Measles, HBV and HCV.



**CDC Guideline & D-125**

D-125 meets CDC selection criteria for Environmental surface disinfection of Health Care facilities

A product of :

**Microgen Hygiene Pvt. Ltd.**

614, Midas, Sahar Plaza, Andheri-Kurla Road, Andheri ( East ). Mumbai - 400 059.

**Tel :** +91 22 2825 1870 **Fax :** +91 22 2832 0335

**Website :** www.microgenindia.co

## Areas of Application :

### Areas of Usage :

Hard, non-porous surfaces in : Hospitals (I.C.U / N.I.C.U / A.K.D / C.C.U / O.T.), Radiology, Labour Department, Isolation Wards, All Medical Devices, Ultrasound Scan Centers, I.V.F. Centers, Nursing homes, Cath labs, Dialysis Clinics, Burn Units, Dental Offices, Laboratory & Blood Banks, Clinics, Home Health Cares, Hospices, Ambulance etc.

### Packing :

Bucket : 200 Wipes  
(wipe size: 15 x 25 cm)

**PRECAUTIONS :** CLOSE THE LID TIGHTLY AFTER EACH USE TO AVOID DRYING